BAS

V.1.5
06/2025

ES | FR | EN | IT | PT | DE

BALANZA ELECTRÓNICA

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
ELECTRONIC BALANCE

BILANCIA ELETTRONICA

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marca propiedad de | est une marque de | trade mark propriety of
marchio di proprietà di | marca de propiedad de | Warenzeichen Proprietät von:

Pol. Empordà Internacional C/ Molló, 3
17469 VILAMALLA - (Girona) SPAIN
T. (34) 972 527 212

GIROPES

ES

INDICE

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EN INDEX

1. PRECAUTIONS	22
2. INTRODUCTION	23
2.1 GENERAL	23
2.2 OVERALL VIEW	23
3. INSTALLATION	24
3.1 UNPACKING	24
3.2 INSTALLATION	24
3.3 LEVELING	24
3.4 POWER CONNECTION	24
3.5 RECHARGEABLE BATTERY OPERATION	25
3.5.1 CHARGING THE BATTERY	25
3.5.2 BATTERY MAINTENANCE	25
3.6 RS232 (9PIN D TYPE CONNECTOR)	25
3.6.1 CONTINUOUSLY OUTPUT PROTOCOL	25
3.6.2 ASK COMMAND	25
4. DESCRIPTION	26
5. OPERATION	27
5.1 POWER ON/OFF	27
5.2 ZERO	27
5.3 TARE	27
5.4 BY10	27
5.5 PERCENT WEIGHING	28
5.6 PARTS COUNTING	28
5.7 ACCUMULATION	28
5.7.1 MEMORY RECALL	29
5.7.2 MEMORY CLEAR	29
5.8 PRINTING SETTING	29
5.9 BACKLIGHT SETTING	29
5.10 RS232 OUTPUT	30
5.11 USB OUTPUT	30
6. PARAMETERS	31
7. CALIBRATION	33
7.1 LINEAR CALIBRATION	33
7.2 NORMAL CALIBRATION	33
8. ERROR CODES AND TROUBLESHOOTING	34
8.1 ERROR CODES	34
8.2 TROUBLESHOOTING	34
9. WARRANTY	35

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1. PRECAUTIONS

	 WARNING
	<p>DISCONNECT ALL POWER TO THIS UNIT BEFORE INSTALLING, CLEANING, OR SERVICING. FAILURE TO DO SO COULD RESULT IN BODILY HARM OR DAMAGE THE UNIT.</p>

 CAUTION
<ul style="list-style-type: none"> - PERMIT ONLY QUALIFIED PERSONS TO SERVICE THE INSTRUMENT. - BEFORE CONNECTING OR DISCONNECTING ANY COMPONENTS, REMOVE THE POWER. - FAILURE TO OBSERVE THESE PRECAUTIONS BODILY HARM OR DAMAGE TO OR DESTRUCTION OF THE EQUIPMENT.

- The weighing scale is a precision electronic instrument, handle it carefully.
- Do not install the scale in direct sunlight.
- Verify the local voltage and receptacle type are correct for the scale.
- Only use original adaptor, other could cause damage to the scale.
- Pluggable equipment must be installed near an easily accessible socket outlet.
- Avoid unstable power sources. Do not use near large users of electricity such as welding equipment or large motors.
- Avoid sudden temperature changes, vibration, wind and water.
- Avoid heavy RF noise.
- Keep the scale clean.

2. INTRODUCTION

2.1 GENERAL

- The BAS series balances, that amplifies signals from a load cell, converts it to digital data and displays it as a mass value.
- It is accurate, fast and versatile series of general purpose balances with counting, % weighing functions and accumulation.
- 16.5 mm LCD with white LED back light display.
- All the keypads are light touch switches.
- Battery provide up to 35 hours of continues use (with out backlight).
- Capacity from 600g to 6000g.
- RS232 port and USB type B port for communication with other devices.
- All units include automatic zero tracking, and an accumulation facility that allows the count to be stored and recalled as an accumulated total.

2.2 OVERALL VIEW

3. INSTALLATION

3.1 UNPACKING

When you receive the scale, inspect it to make sure that it is not damaged and that all parts are included.

- Remove the scale from the carton.
- Remove the protective covering. Store the packaging to use if you need to transport the scale later.
- Inspect the scale for damage.
- Make sure all components are included:
 1. Scale.
 2. Adaptor.
 3. Manual.
 4. Stainless steel pan.
 5. Wind shield top cover (specific model).

3.2 INSTALLATION

Remove transport lock from bottom of the scale base.

Install weighing pan as follows:

Models **BAS 600 / 1200 / 3000**


1. Place the scale on a level surface.
2. Install the plastic pan into holder and fix the pan with screw.
3. Install the stainless steel pan on top of the plastic pan.

Models **BAS 6000**

1. Place the scale on a level surface.
2. Install the plastic pan into holder and fix the pan with screw.
3. Install the stainless steel pan on top of the plastic pan.

3.3 LEVELING

- The scale is equipped with a bubble level located left side of weight display.
- Use the adjustable leveling feet located on the bottom of the scale until the bubble appears in the center of the indicator.

3.4 POWER CONNECTION

- Connect the adaptor pin in to the adaptor jack.
- Adaptor jack is locating beside side housing of the scale.
- Adaptor connects into your AC power socket. Pluggable equipment must be installed near an easily accessible socket outlet with a protective ground/ earth contact.

3.5 RECHARGEABLE BATTERY OPERATION

NOTE: Please charge the battery before using the scale for the first time.

- A symbol on the weight display indicates that the battery is getting low, it's time to charge the battery with the AC power. Approximately 1 hour of instrument usage is left; afterwards it will shut off automatically.
- Please use the supplied battery charger for charging the battery immediately, or scale cannot be used.


3.5.1 CHARGING THE BATTERY

- Before the first use, the battery should be charged by connecting it to the mains power supply for at least 1.5 hours.
- When the scale is plugged into the mains power the internal battery will be recharged.
- Verify that the AC power socket outlet is properly protected.

3.5.2 BATTERY MAINTENANCE

- Do not use any other type of power adaptor than the one supplied with the scale.
- Recharge battery every three months when not in use.
- If the scale is not used for an extended period of time, remove the battery from the battery compartment to avoid leakage.
- Store the battery in a sealed bag or box in a dry, temperate environment.

3.6 RS232 (9PIN D TYPE CONNECTOR)

Pin 2	RXD	Input	Receiving data
Pin 3	TXD	Output	Transmission data
Pin 5	GND	-	Ground signal

9pin D Connector

Indicator	Computer / Printer
Pin 2	Pin 3
Pin 3	Pin 2
Pin 5	Pin 5

NOTE: If data is not printed, inter-change the pin 2 and pin 3 connections.

3.6.1 CONTINUOUSLY OUTPUT PROTOCOL

HEADER1: ST=STABLE, US=UNSTABLE

HEADER2: NT=NET, GS=GROSS

WEIGHING UNIT : g

Data Communication example

S	T	,	G	S	,	-/ □								g			CR	LF
-HEADER-			-HEADER 2-				-WEIGHT DATA-					-UNIT-		TERMINATOR				

3.6.2 ASK COMMAND

Command	Function	Pinout examples
W	Weighing value for the weight (stable or unstable) is sent via the RS232 interface	US,GS, 588g ST,GS, 588g
T	No data are sent, the balance carries out the tare function.	-
Z	No data are sent, the zero-display appears.	-

4. DESCRIPTION

KEYBOARD

- Print the results.
- Exit menu or back to weighing mode.

On/Off key, used to turn the scale on or off.

- Mode key, used for extended or auxiliary mode depending of the scale capacity model.
- In setup menu, used to scroll forward the menu and sub-menu.

- Count key, used to enter the counting operation.
- In setup menu, used to scroll downward the menu and sub-menu.

Percentage key, used to set percentage weighing function.

Tare key, used to perform a tare function, Subtracts weights.

- Zero key, used to clear and zero the display.
- Enter key, used to enter the selected menu, sub-menu and setting.

DISPLAY	FUNCTION	DISPLAY	FUNCTION
→0←	Indicator for Zero display	mom Kg ㄗㄗㄗ	Indicator for units
→T←	Indicator for Tare display	%	Indicator for Percent weighing
00000000000000000000	Indicator for weighing capacity graph	[Battery icon]	Indicator for charging status of battery
NET	Indicator for Net weight	[Battery icon]	Voltage has dropped
◎	Indicator for Display stability	[Battery icon]	Low Voltage
Pcs	Indicator for piece counting	[Battery icon]	Fully Charged

5. OPERATION

Initial Start-up

Warm-up time of 30 minutes stabilizes the measured values after switching on.

5.1 POWER ON/OFF

Switch on the scale by pressing . The display will show "warm up" prompt for 40 seconds with the dynamic graphic bar, the motor noise of the loading system for the internal adjustment weight can be heard. The balance carries out a self-test, as soon as the weight display appears, the balance is ready for weighing.

If you want to switch off press the key again.

5.2 ZERO

Environmental conditions can cause the scale to not reach exactly zero even if the pan is not subjected to any stress. You can set the display of your balance to zero any time by pressing  key and ensure that the weighing starts at zero.

5.3 TARE

The weight of any container can be tared by pressing  button so that with subsequent weighing the net weight of the object being weighed is always displayed.

- Load weight on the pan.
- Press key. Zero is displayed, and tare is subtracted.
- Remove weight from the platform. Tared weight is displayed. Press  key. Zero is displayed, tare weight is cleared.

5.4 BY10

Press  key in weighing mode to toggle the display view between the normal display image, the auxiliary indicating device, or the extended indication display, this depends on the model by the scale's maximum capacity.

5.5 PERCENT WEIGHING

The scale can set a sample weight to be shown as 100%. Then any other weights place on the scale, it will be displayed as a percentage of the original sample.

For example:

1. 350g weight place on the scale and follow by press the  key is pressed the display will show 100.00%.
2. Remove the weight, and ensure display is zero.
3. Place 300g weight on the platform, display will show 85.71 % as for the percentage of 350g (100%).
4. The weighing may be amended on the basis of greater numbers of samples, improves the accuracy of percentage large quantities.
5. Press  key again to back to normal weighing mode.

5.6 PARTS COUNTING

Press  key enter the parts counting mode and select the counting options by pressing . Display will show:

(FIRST, LOAD THE SAMPLE WEIGHT)

sp 10	For 10 pieces
sP 20	For 20 pieces
sP 50	For 50 pieces
sP 100	For 100 pieces
sP 200	For 200 pieces

Select the counting option and press  to confirm. Then can add more weight display will be the weight show the number of parts.

Parts Counting Operation

If necessary place a container on the platform and press  to make zero.

1. Place the load on the platform.
2. When stable indicator on, press  key enter the parts counting mode.
3. Press  key to the parts quantity as per the option.
4. Press  key to confirm, display will show the quantity.
5. Then can add goods on the platform, display will update the parts quantity automatically.
6. Press  key to change to normal weighing mode, when in counting mode.

5.7 ACCUMULATION

The scale can be set to accumulate manually by pressing  key.

Before operation scale should be stable and return to zero, accumulation available only when weight more than 20d. To enable accumulation function, select parameter:

F2 COM > RS232 > PRINT > TPuP or LP50

F3 ACC > ON

NOTE: CAL switch needs to be pressed.

Accumulation Operation

- Place the load on the platform.
- Press  key, when displayed STABLE indication.
- Display will be show [ACC !]. Then will be show the total saved value. These displays will be shown only three seconds. If the RS-232 interface is installed the weight data will be send to printer.
- Remove the weight from the pan.
- When display get zero and stable then place the second weight.
- It can continue until the memory gets fully.

5.7.1 MEMORY RECALL

To recall the memory press  key in zero point.

Display will show [ACC X].

(X: Total number of accumulation) then will show the total saved value. These displays will be shown only three seconds.

5.7.2. MEMORY CLEAR

To clear the memory, press  key to check and press  during the accumulation displayed. Display will show weighing display. All accumulation memory cleared from the memory.

5.8 PRINTING SETTING

To enable the printer function, see parameters:

F2 COM > RS232 > BAUD RATE > 9600

F2 COM > RS232 > PRINT > TPuP or LP50

TPUP is for receipt printer, and LP50 is for label printer:

1. Turn on the scale.
2. Press  key during self checking, display will show **F1 bl**.
3. Press  key until display show **F2 COM**.
4. Press  key to confirm, display show **RS 232**.
5. Press  key to enter, display show **PRINT**.
6. Press  key to select **LP50**. If necessary press  key to select **TPuP** printer.
7. Press  key to confirm, display show **PRINT**.
8. Press  key to select **BAUD RATE**. Press  key to confirm.
9. Press key  to select **b 9600** and use  key to confirm.
10. Press  key 3 times to back to normal weighing mode.
11. Load the weight on the platter, after stable indicator on.
12. Press  key. The format will be sent to print out.

5.9 BACKLIGHT SETTING

1. Press  key during self checking to select parameter **F1 bl**.
2. Press  key to select (**BL AU/ BL ON/ BL OFF**).
3. Press  key to confirm.



F1bl	bl Au	To set auto option. When start to use, back light will turn on and when stop the operation, back light will turn off.
	bl on	To set always on. After turn on the power, back light will be on.
	bl off	To set back light turn off. No back light in the operations.

NOTE: When battery gets low, the backlight function will not be available.

5.10 RS232 OUTPUT

To enable the communication of balance, it is required to enter parameters to select: **F2 COM > RS232 > PROTCL > CONT.**

Operation:

1. Turn on the scale.
2. Press  key during self checking, display will show **F1 bL**.
3. Press  key until display show **F2 COM**.
4. Press  key to confirm, display show **RS 232**.
5. Press  key to enter, display show **PRINT**.
6. Press  key to select **PROTCL**.
7. Press  key to enter, display show **OFF**.
8. Press  key until display show **CONT**.
9. Press  key to confirm, display show **PROTCL**.
10. Press  key to select **BAUD RATE** and use  key to confirm.
11. Press  key to select **b 9600** and use  key to confirm.
12. Press  key 3 times to back to normal weighing mode.
13. Load the weight on the platter, the weight value will be sent to the computer.

5.11 USB OUTPUT

The balance is equipped with an USB port (type B) for data output.

Preparing a USB cable to connect with PC or printer, one side of cable connect with computer, the other side of cable insert to the USB port which is located on the back side of housing.



Scan the QR code below to download the USB driver and view USB installation and operation guide:

USB Driver

To enable the communication of balance, it's required to enter parameters to select **F2 COM > USB > CONT.**

Operation:

1. Turn on the scale.
2. Press  key during self checking, display will show **F1 bL**.
3. Press  key until display show **F2 COM**.
4. Press  key to confirm, display show **RS 232**.
5. Press  key to select **USB**.
6. Press  key to enter, display show **PRINT**.
7. Press  key to select **PROTCL**.
8. Press  key to enter, display show **OFF**.
9. Press  key until display show **CONT**.
10. Press  key to confirm, display show **PROTCL**.
11. Press  key to select **BAUD RATE** and use  key to confirm.
12. Press  key to select **b 9600** and use  key to confirm.
13. Press  key 3 times to back to normal weighing mode.
14. Load the weight on the platter, the weight value will be sent to the computer.

6. PARAMETERS

Enter the Menu

- Turn on the scale. Press  during that start up until display will be shown **F1 bL**.

Choose the Menu

- Press  or , it can choose menu block or options one by one.

Enter the Selected Menu

- Press , it can confirm which will be shown displayed.

Enter in to TECH

- When display showed **F1 r n**, press ,  and  keys to enter the function.

Escape from the Menu

- Press  key, it can escape from the menu to weighing mode.


MENU	SUBMENU	DESCRIPTION	DEFAULT		
F1 bL	bL Au	To set back light automatic on	BL AU		
	bL on	To set back light always on			
	bL off	To set back light always off			
F2 con	r5232	Pr int	oFF	Printer function off	OFF
		LP50	To connect LP50 label printer		
		tPUP	To connect Tpup receipt printer		
	Protcl	cont	Send data continuous	CONT	CONT
		ASK	ASK mode Command W: Read weight Command T: Tare Command Z: Zero		
		KcP	KCP protocol		
		oFF	Communication port function off		
	bAud	4800	To set baud rate	9600	
		9600			
		57600			
		115200			
	usb	Pr int	oFF	Printer function off	OFF
			LP50	To connect LP50 label printer	
tPUP			To connect Tpup receipt printer		

	usb	Protcl	cont	Send data continuous	
			ASK	ASK mode Command W: Read weight Command T: Tare Command Z: Zero	
			KCP	KCP protocol	
			oFF	Communication port function off	
	bAud	4800	To set baud rate	9600	
		9600			
		57600			
		115200			
	LAB	LAB 1	Format 1		
		LAB2	Format 2		
F3Acc	Acc on		Accumulation function enable	ON	
	Acc oFF		Accumulation function disable		
F4n-t	Sbt ich	on	Key delay function enable	ON	
		oFF	Key delay function disable		
	t iNE	3	Key delay function is suitable for unstable environments. Select the key delay time (3 / 5 / 8 seconds)	3	
		5			
8					
F5btt	on		To enable optimized weighing experience	OFF	
	oFF		To disable optimized weighing experience		
tEch	P in		Enter the password		
P 1L in			Linear Calibration		
P2cAL			Normal Calibration		
P3cnt	HHHH		This display will show XXXXX for indicating the internal counts		
P4GrA	HHHHH		To set local gravity Press % key to move digits right,press TARE key to increment the digits		
P5cAP	150G/300G/600G/ 1500G/ 3000G/6000G		To set Capacity		
P8SPd	LoB		Set the ADC speed	LOW	
	N id				
	h iGh				



7. CALIBRATION

Before entering the tech menu, it is required to press the calibration switch, which is located below the scale.

7.1 LINEAR CALIBRATION

- Turn on the scale.
- Press  key during self checking, display will show $F \quad IbL$.
- Press  key until display will show $tEcH$.
- Press  key to enter calibration, display will show $P \quad r$.
- In succession press ,  and  keys to enter the function, display will show $P \quad IL \quad r$.
- Press  key to enter calibration, display will show $P \quad r$.
- In succession press , ,  keys to enter the function display will shown $uRLoAd$.
- Ensure the platform is empty.
- 3 seconds later, display will show 1/6 of the capacity.
- Place the first calibration (1/6 of the capacity) mass weight on the pan.
- 3 seconds later, display will show $uRLoAd$.
- Remove the first calibration (1/6 of the capacity) mass weight from the pan.
- 3 seconds later, display will show 2/6 of the capacity.
- Place the second calibration (2/6 of the capacity) mass weight on the pan.
- 3 seconds later, display will show $uRLoAd$.
- Remove the second calibration (2/6 of the capacity) mass weight from the pan.
- 3 seconds later, display will show 3/6 of the capacity
- Place the third calibration (3/6 of the capacity) mass weight on the pan.
- 3 seconds later, display will show $uRLoAd$
- Remove the third calibration (3/6 of the capacity) mass weight from the pan.
- 3 seconds later, display will show 4/6 of the capacity.
- Place the fourth calibration (4/6 of the capacity) mass weight on the pan.
- 3 seconds later, display will show $uRLoAd$.
- Remove the fourth calibration (4/6 of the capacity) mass weight from the pan.
- 3 seconds later, display will show 5/6 of the capacity.
- Place the fifth calibration (5/6 of the capacity) mass weight on the pan.
- 3 seconds later, display will show $uRLoAd$.
- Remove the fifth calibration (5/6 of the capacity) mass weight from the pan.
- 3 seconds later, display will show full capacity.
- Place the sixth calibration (full capacity) mass weight on the pan.
- 3 seconds later, display will show $uRLoAd$.
- Remove the sixth calibration (full capacity) mass weight from the pan.
- 10 seconds later, scale will show $-PAsS-$ and back to setting mode.

7.2 NORMAL CALIBRATION

- Turn on the scale.
- Press  key during self checking, display will show $F \quad IbL$.
- Press  key until display will show $tEcH$.
- Press  key to enter calibration, display will show $P \quad r$.
- In succession press ,  and  keys to enter the function, display will show $P \quad IL \quad r$.
- Press  key to select $P \quad cAL$.
- Press  key to enter calibration, display will show $uRLoAd$.
- Remove all the weight from the platform.
- 3 seconds later, display will show the mass weight.
- If necessary press  key to select the weight as you desired.
- Place the calibration weight on the platform, a few seconds later, scale will show $-PAsS-$ and back to setting mode.

8. ERROR CODES AND TROUBLESHOOTING

8.1 ERROR CODES

Error Message	Description	Solution
--oL--	Maximum load exceed	Descargar o reducir el peso
Err 4	Zero setting error	*Zero setting range exceeded due to pressing ZERO key. (20%max) *Initial zero setting range exceeded due to switching on balance (4%max). Make sure platform empty.
Err 7	Sample or Percentage Error	Please check the load weight, must be more than 20d
--uL--	Underload	Minus weight, check the platform and restart or calibrate.
FA iL	Calibration Error	Check the test weights & Re calibrate

8.2 TROUBLESHOOTING

Problem	Possible cause	Common solutions
Display is blank, No turn on test	Mains power is turned off Power supply not plugged in Power supply faulty Internal battery not charged Display turned off	Check power is going to the scale and switch is on. Verify the voltage going to the scale matches the power supply labels on the power module or scale.
Display blank after turn on test, Error message Display is locked	Pan not installed Unstable weight Load cell damaged Mechanics damaged	Check the pans are installed correctly. Try turning the scale on again.
OL or FULL appears on display	Maximum capacity exceeded Load Cell or mechanics damaged Power supply faulty	Check the pans are installed correctly. Try turning the scale on again.
Display is unstable	Drafts or air currents Obstruction under pan Sample is moving (animal weighing) Vibrations through table Temperature changed dramatically Power supply faulty	Verify the scale is in a acceptable location and on a good table. Verify the power supply is correct for the scale.
Weight value incorrect	Calibration error, Recalibrate Linearity error, set Linearity Unit calibrated with inaccurate weight Balance not level Obstruction between sample and cover Wrong unit of weight displayed	Calibrate again, paying special attention to the mass used, the stability of the scale, and the weighing units required. If linearity can be set by the user it will be described in the user manual. Check pan is installed correctly. Verify installation is acceptable. The display of the balance is not at zero. Warm-up time was ignored.
Cannot use Full Capacity	Over load Stops hitting pan support or hitting bottom of load cell Shipping screw not removed if applicable Electronic problem on A/D Parameters set incorrectly Load Cell or mechanics Damaged	Look for obstruction under pan, shipping screws, and verify pan installation. Check the weighing units used.

Not Linear	Overload stops hitting too soon Load cell or mechanics damaged A/D damaged	Look for obstruction under pan, shipping screws, and verify pan installation. If linearity can be set by the user it will be described in the user manual.
Off Center Loading error	Adjust mechanics Overload Stops not correct Load Cell damaged	Look for obstruction under pan, shipping screws, and verify pan installation.
Battery will not charge	Mains voltage not present or too low Charging circuit failure Battery Failure	Verify the batteries are rechargeable types. Check power supply voltage is correct

9. WARRANTY

GIROPES S.L. with its mark BAXTRAN guarantees its products, bought through an authorized channel, with reference to manufacturing or material defects during a period of two years from the delivery date, according to the current legislation, except when there is a specific agreement or in the acceptance of the order.

- » The warranty does not include damages produced by the wear of the material due to its use, inappropriate conservation or maintenance, wrong storage or use, modifications introduced without the consent, on paper, from Giropes, and in general due to causes beyond the control of our company, Giropes.
- » The reparations will be done by Giropes. The customer must pay for the dismantling, the packaging, the transport, the taxes, custom duty... produced by the delivery of the material to Giropes and its delivery to the customer. The customer must do the delivery of the material via the RMA system which Giropes has.
- » Giropes can agree with the customer the reparation or the replacement of the faulty pieces. In such case, the expenses of the journey of the members of the technical service will be paid by the customer.
- » In case that it is needed to replace a faulty product completely, all the expenses will be paid by BAXTRAN, including the arriving and shipping dues.
- » Giropes is not going to pay for the reparations done by third parties.
- » The repairing or replacement of a faulty piece does not vary the guarantee period of the supplied product. However, the piece that has been repaired or replaced has a guarantee period of two years from the date when the reparation or replacement has been done.
- » Giropes is not responsible, in any case, of indirect damages and/or as a consequence of the supplying.
- » The Customer will be responsible of getting rid of the elements supplied by BAXTRAN, following the necessary regulation rules, by its own account and responsibility.
- » With regard to product consumables (battery, adapter, thermal heads, etc.) the warranty will be as follows:
 - » Free change of the consumable during the first 6 months from the moment of delivery of the product.
 - » Free replacement of consumables from 6 to 24 months from delivery of the product provided that the product is proven to have been used properly prior to the failure of the product.

WARRANTY FOR DAMAGES DONE DURING THE TRANSPORT:

- » Customers out of Spain but in European Union have 7 days from the date of delivery of the merchandise to communicate to BAXTRAN any type of damage produced by the transport in case the transport has been hired by Baxtran. If the complaining is not done during this period, BAXTRAN will not be responsible and in charge to pay for any type of repairing or substitution if the product has been damaged. This period could change by other country's law.
- » We remind you that there is at your disposal a section called "Reservations entry" to specify, when signing the delivery note that the package contains damages such as hits, breaks or others that can affect the product and cannot be seen without opening the package. We strongly recommend you to complete with your observations the delivery note from the deliverer before signing it. All the members of our staff have an EC Certify and the user's guide for our customers.

RMA

GIROPES announced the introduction of a RMA service (Return Merchandise Authorization) with the objective to improve steps for products returns and/ or repair and give a better after-sales services. From now, RMA number is going to be compulsory to start any Technical Service.

If you are not familiar with our RMA system, please click the link and follow the instructions: <https://www.baxtran.com/en/rma.html>

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